

Computer Setup

Set Trusted Sites Custom Level

Setting the trusted sites custom level enables security settings so that the AdvancedMD software can be installed at your workstation.

- 1. Launch Internet Explorer 32-bit (latest 2 versions are compliant)
- 2. Go to www.onpointmedicalsolutions.com
- 3. From the **Tools** menu, select **Internet Options**.
- 4. Select the **Security** tab.
- 5. Click Trusted sites.
- 6. Click the Sites button.
- 7. Type https://*.advancedmd.com and click Add.
- 8. Uncheck the box, "Requires server verification (https)."
- 9. Type http://*.advancedmd.com and click Add.
- 10. Type http://*.onpointmedicalsolutions.com and click Add.
- 11. Click **OK**.
- 12. Click the **Custom Level** button.
- 13. Reset custom settings.
 - a. Reset to LOW
 - b. Click Reset
- 14. Click **OK**.
- 15. Click **OK**.

Allow Pop ups

This system uses pop up screens. Pop up blocker needs to be set up to allow pop ups from the AdvancedMD system. In addition to these Internet Explorer settings any third party toolbars such as Yahoo, Google, Bing etc. must also be removed or the pop-up blockers disabled.

- 1. From the **Tools** menu, point to **Pop-up Blocker**, and then select **Pop-up Blocker Settings**. Pop-up blocker setting may also be in the **Privacy tab**.
- 2. Type *.advancedmd.com and click Add.
- 3. In the **Blocking Level** box select **Low**.
- 4. Click Close
- 5. Close all instances of Internet Explorer This is necessary for settings to take effect.



Add Onpoint - Log In Shortcut ICON on Desktop

- 1. Launch Internet Explorer 32-bit (latest 2 versions are compliant)
- 2. Go to www.Onpointmedicalsolutions.com
- 3. Go to CLIENT LOGIN
- 4. Go to Favorites
 - a. Click Add to Favorites
 - b. Rename favorite to **Onpoint Log In**
 - c. Click Add
- 5. Open Favorites
 - a. Right Click on Onpoint Log In, select copy
 - b. Right click on your desktop, select paste shortcut

Log in for the first time

- 1. Launch **Onpoint Log In** shortcut
- 2. **Enter username, password and office key** This information is provided to you by your Onpoint Trainer or representative and is necessary to complete the install process.
- 3. Click Onpoint PM Here you will be launched in to the Onpoint / AdvancedMD software.
- 4. Click **Utilities**
- 5. Select ActiveX controls
- 6. Select Test and Install controls
- Your results should indicate success. If not contact your Onpoint Trainer or representative.